

ALKBH5 Antibody

Catalog No: #44852

Package Size: #44852-1 50ul #44852-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ALKBH5 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	ALKBH5antibody detects endogenous levels of total?ALKBH5
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ALKBH5
Target Name	ALKBH5
Other Names	ABH5 antibody AlkB, alkylation repair homolog 5 (E. coli) antibody AlkB, alkylation repair homolog 5 antibody Alkylated DNA repair protein alkB homolog 5 antibody Alpha ketoglutarate dependent dioxygenase alkB homolog 5 antibody OFOXD antibody OFOXD1 antibody Oxoglutarate and iron-dependent oxygenas
Accession No.	Swiss-Prot#: Q6P6C2NCBI Gene ID:54890
Uniprot	Q6P6C2
GeneID	54890;
Calculated MW	44 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C


Application Details

Western blotting: □ 1:500-1:2000

Images



Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).

Note: This product is for in vitro research use only