

BCL7B Antibody

Catalog No: #44853

Package Size: #44853-1 50ul #44853-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	BCL7B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	BCL7B antibody detects endogenous levels of total BCL7B
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human BCL7B
Target Name	BCL7B
Other Names	B cell CLL/lymphoma 7 protein family member B antibody B cell CLL/lymphoma 7B antibody B-cell CLL/lymphoma 7 protein family member B antibody BCL 7B antibody BCL7B antibody BCL7B_HUMAN antibody
Accession No.	Swiss-Prot#: Q9BQE9 NCBI Gene ID:9275
Uniprot	Q9BQE9
GeneID	9275;
Calculated MW	22 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500-1:2000

Images



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only