

## PCTK3 Antibody

Catalog No: #44887

Package Size: #44887-1 50ul #44887-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	PCTK3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Immunogen Type	Peptide
Immunogen Description	Fusion protein of human PCTK3
Target Name	PCTK3
Other Names	Cdk18 antibody CDK18_HUMAN antibody Cell division protein kinase 18 antibody  Cyclin-dependent kinase 18 antibody PCTAIRE 3 antibody PCTAIRE antibody PCTAIRE motif protein kinase 3 antibody PCTAIRE protein kinase 3 antibody PCTAIRE-motif protein kinase 3 antibody
Accession No.	Swiss-Prot#: Q07002NCBI Gene ID:5129
Uniprot	Q07002
GeneID	5129;
Calculated MW	54 kDa
Concentration	1.0mg/mL
Formulation	pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
Storage	Store at -20°C

## Application Details

Western blotting: □ 1:500-1:2000

## Images



Gel: 8%SDS-PAGE, Lysate: 40  $\mu$ g, Lane: HL-60 cell lysate,  
Primary antibody: PCK3 Antibody at dilution 1/250,  
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution,  
Exposure time: 5 minutes

---

Note: This product is for in vitro research use only