

CASP6 Conjugated Antibody

Catalog No: #C32434



Package Size: #C32434-AF350 100ul #C32434-AF405 100ul #C32434-AF488 100ul
 #C32434-AF555 100ul #C32434-AF594 100ul #C32434-AF647 100ul
 #C32434-AF680 100ul #C32434-AF750 100ul #C32434-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	CASP6 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous level of total CASP6 protein.
Immunogen Description	Recombinant protein of human CASP6.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	MCH2;CASP6;Caspase-6
Accession No.	Swiss-Prot#:P55212NCBI Gene ID:839
Uniprot	P55212
GeneID	839;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	22
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Product Description

Antibodies were purified by affinity purification using immunogen.

Background

Caspase-6 (Mch2) is one of the major executioner caspases functioning in cellular apoptotic processes (1,2). Upon apoptotic stimulation, initiator caspases such as caspase-9 are cleaved and activated (3). The activated upstream caspases further process downstream executioner caspases, such as caspase-3 and caspase-6, by cleaving them into large and small subunits, thereby initiating a caspase cascade leading to apoptosis (4,5). One of the major targets for caspase-6 is the membrane associated protein lamin A (6). The cleavage of this protein causes cell membrane malfunction, membrane blebbing and eventual cell death.

Note: This product is for in vitro research use only