

## TAOK2 Antibody

Catalog No: #44892

Package Size: #44892-1 50ul #44892-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	TAOK2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	TAOK2antibody detects endogenous levels of total?TAOK2
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human TAOK2
Target Name	TAOK2
Other Names	1110033K02Rik antibody B230344N16 antibody hKFC C antibody  hKFC-C antibody KIAA0881 antibody Kinase from chicken homolog C antibody MAP3K17 antibody mKIAA0881 antibody Prostate derived STE20 like kinase 1 antibody Prostate derived STE20 like kinase PSK antib
Accession No.	Swiss-Prot#: Q9UL54NCBI Gene ID:9344
Uniprot	Q9UL54
GeneID	9344;
Calculated MW	120 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: □ 1:500-1:2000

## Images



Western Blot analysis of HEK293 lysates, using primary antibody at 1:1000 dilution.

---

Note: This product is for in vitro research use only