

PRDX6 Conjugated Antibody

Catalog No: #C32561



Package Size: #C32561-AF350 100ul #C32561-AF405 100ul #C32561-AF488 100ul
 #C32561-AF555 100ul #C32561-AF594 100ul #C32561-AF647 100ul
 #C32561-AF680 100ul #C32561-AF750 100ul #C32561-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PRDX6 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous level of total PRDX6 protein.
Immunogen Description	Recombinant protein of human PRDX6.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	1-Cys;AOP2;KIAA0106;MGC46173;NSGPx
Accession No.	Swiss-Prot#:P30041NCBI Gene ID:9588
Uniprot	P30041
GeneID	9588;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	25
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Product Description

Antibodies were purified by affinity purification using immunogen.

Background

Prdx6(Peroxiredoxin-6) , a unique member of the peroxiredoxin family of antioxidants, is highly expressed in liver and protects cells from oxidative damage by reducing H₂O₂ and various lipid Peroxides(PMID: 17382207).It is also named as AOP2, KIAA0106. It can form a dimer(PMID:20500660).PRDX6 is expressed in all major organs, with a particularly high level in lung(PMID:15890616).This antibody is specific to PRDX6.

Note: This product is for in vitro research use only