

GPR12 Antibody

Catalog No: #44910



Package Size: #44910-1 50ul #44910-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GPR12 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	GPR12antibody detects endogenous levels of total?GPR12
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human GPR12
Target Name	GPR12
Other Names	G protein coupled receptor 123 antibody GP123_HUMAN antibody GPR123 antibody KIAA1828 antibody Probable G-protein coupled receptor 123 antibody
Accession No.	Swiss-Prot#: P47775NCBI Gene ID:2835
Uniprot	P47775
GeneID	2835;
Calculated MW	50 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500-1:2000

Note: This product is for in vitro research use only