

GPR37 Antibody

Catalog No: #44926

Package Size: #44926-1 50ul #44926-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GPR37 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	GPR37antibody detects endogenous levels of total?GPR37
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human GPR37
Target Name	GPR37
Other Names	EDNRBL antibody Endothelin B receptor like protein 1 antibody Endothelin B receptor-like protein 1 antibody Endothelin receptor type B like antibody ETBR LP 1 antibody ETBR-LP-1 antibody ETBRLP1 antibody G protein coupled receptor 37 antibody G protein coupled re
Accession No.	Swiss-Prot#: O15354NCBI Gene ID:2861
Uniprot	O15354
GeneID	2861;
Calculated MW	67 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C


Application Details

Western blotting: □ 1:500-1:2000

Images



Western blot analysis of GPR37 Antibody. The lane on the right is blocked with the GPR37 peptide.



Western blot analysis of the lysates from HT-29 cells using GPR37 antibody.

Note: This product is for in vitro research use only