

GPR61 Antibody

Catalog No: #44933

Package Size: #44933-1 50ul #44933-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	GPR61 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	GPR61 antibody detects endogenous levels of total GPR61
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human GPR61
Target Name	GPR61
Other Names	BALGR antibody Biogenic amine receptor like G protein coupled receptor antibody Biogenic amine receptor-like G-protein coupled receptor antibody G protein coupled receptor 61 antibody GPCR3 antibody Gpr61 antibody GPR61_HUMAN antibody Probable G-protein coupled recep
Accession No.	Swiss-Prot#: Q9BZJ8 NCBI Gene ID:83873
Uniprot	Q9BZJ8
GeneID	83873;
Calculated MW	49 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500-1:2000

Note: This product is for in vitro research use only