

## GPR81 Antibody

Catalog No: #44942

Package Size: #44942-1 50ul #44942-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	GPR81 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	GPR81 antibody detects endogenous levels of total GPR81
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human GPR81
Target Name	GPR81
Other Names	FKSG80 antibody FKSG80 protein antibody G protein coupled receptor 104 antibody  G protein coupled receptor 81 antibody G protein coupled receptor GPR81 antibody G-protein coupled receptor 104 antibody G-protein coupled receptor 81 antibody GPR 104 antibody GPR 8
Accession No.	Swiss-Prot#: Q9BXC0 NCBI Gene ID:27198
Uniprot	Q9BXC0
GeneID	27198;
Calculated MW	39 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: □ 1:500-1:2000

## Images



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only