

## GPR89A Antibody

Catalog No: #44947

Package Size: #44947-1 50ul #44947-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	GPR89A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	GPR89Antibody detects endogenous levels of total?GPR89A
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human GPR89A
Target Name	GPR89A
Other Names	G protein coupled receptor 89A antibody Golgi pH regulator A antibody Golgi pH regulator B antibody  GPHR antibody GPR89 antibody GPR89B antibody Protein GPR89 antibody Putative MAPK activating protein PM01 antibody Putative NF kappa B activating protein 90 antib
Accession No.	Swiss-Prot#: Q9Y302NCBI Gene ID:51463/653519
Uniprot	B7ZQA6
GeneID	51463;653519;
Calculated MW	53 kDa
Concentration	1.0mg/mL
Formulation	Reconstitute in 100ul of sterile H2O with 50% glyceroloO 0.01% thiomersal.
Storage	Store at -20°C

## Application Details

Western blotting: □ 1:500-1:2000

Note: This product is for in vitro research use only