

GNRHR2 Antibody

Catalog No: #44960

Package Size: #44960-1 50ul #44960-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GNRHR2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	GNRHR2 antibody detects endogenous levels of total GNRHR2
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human GNRHR2
Target Name	GNRHR2
Other Names	GnRH IIR antibody Putative gonadotropin-releasing hormone II receptor antibody Type II GnRH receptor antibody
Accession No.	Swiss-Prot#: Q96P88 NCBI Gene ID:114814
Uniprot	Q96P88
GeneID	114814;
Calculated MW	42 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500-1:2000

Images



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only