

GPR141 Antibody

Catalog No: #44962



Package Size: #44962-1 50ul #44962-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GPR141 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	GPR141 antibody detects endogenous levels of total GPR141
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human GPR141
Target Name	GPR141
Other Names	GPR141 antibody G protein coupled receptor 141 antibody G-protein coupled receptor PGR13 antibody PGR13 antibody Probable G-protein coupled receptor 141 antibody
Accession No.	Swiss-Prot#: Q7Z602 NCBI Gene ID:353345
Uniprot	Q7Z602
GeneID	353345;
Calculated MW	35 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500-1:2000

Images



Western blot analysis of lysates from HCT116 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only