

LPHN3 Antibody

Catalog No: #44975

Package Size: #44975-1 50ul #44975-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	LPHN3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	LPHN3antibody detects endogenous levels of total?LPHN3
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human LPHN3
Target Name	LPHN3
Other Names	Calcium independent alpha latrotoxin receptor 3 antibody Calcium-independent alpha-latrotoxin receptor 3 antibody CIRL 3 antibody CIRL-3 antibody CIRL3 antibody cl3 latrotoxin receptor antibody KIAA0768 antibody Latrophilin 3 antibody Latrophilin homolog 3 (cow)
Accession No.	Swiss-Prot#: Q9HAR2NCBI Gene ID:23284
Uniprot	Q9HAR2
GeneID	23284;
Calculated MW	162 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500-1:2000

Images



Western Blot analysis of HEK293 lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only