

KIFC1 Antibody

Catalog No: #45025

Package Size: #45025-1 50ul #45025-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	KIFC1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB, IH, IF/IC
Species Reactivity	Hu
Specificity	KIFC1 Antibody detects endogenous levels of total KIFC1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human KIFC1
Target Name	KIFC1
Other Names	KIFC1, HSET, KNSL2, Kinesin-like 2, Kinesin-related protein, Kinesin family member C1, Kinesin-like protein KIFC1, Kinesin-like protein 2, Kinesin-related protein HSET
Accession No.	Swiss-Prot#: Q9BW19NCBI Gene ID:3833
Uniprot	Q9BW19
GeneID	3833;
Calculated MW	74kDa
Concentration	1.0mg/mL
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/2000), IH (1/50 - 1/200), IF/IC (1/50 - 1/200)

Note: This product is for in vitro research use only