

## LRP6 Antibody

Catalog No: #45043

Package Size: #45043-1 50ul #45043-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	LRP6 Antibody
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	LRP6 Antibody detects endogenous levels of total LRP6
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human LRP6
Target Name	LRP6
Other Names	LRP6, ADCAD2, LRP-6
Accession No.	Swiss-Prot#: O75581NCBI Gene ID:4040
Uniprot	O75581
GeneID	4040;
Calculated MW	180 KD
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only