

CAPER Antibody

Catalog No: #45057



Package Size: #45057-1 50ul #45057-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | CAPER Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | Immunogen affinity purified |
| Applications | WB |
| Species Reactivity | Hu Ms Rt |
| Specificity | CAPER Antibody detects endogenous levels of total CAPER |
| Immunogen Type | Peptide |
| Immunogen Description | A synthesized peptide derived from human CAPER |
| Target Name | CAPER |
| Other Names | ATP6AP2, APT6M8-9, ATP6M8-9, ELDF10, HT028, MSTP009, Renin receptor, M8-9, N14F, Renin/prorenin receptor, ATP6IP2, V-ATPase M8.9 subunit, CAPER, MRXE, XMRE |
| Accession No. | Swiss-Prot#: O75787NCBI Gene ID:10159 |
| Uniprot | O75787 |
| GeneID | 10159; |
| Calculated MW | 42 kDa |
| Concentration | 1.0mg/mL |
| Formulation | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | Store at -20°C |

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only