

## Mitofusin Antibody

Catalog No: #45059

Package Size: #45059-1 50ul #45059-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Mitofusin Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	Mitofusin Antibody detects endogenous levels of total Mitofusin
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human Mitofusin
Target Name	Mitofusin
Other Names	MFN1, Fzo homolog, Hfzo2, Mitofusin-1, Hfzo1, Mitofusin 1, Putative transmembrane GTPase, Transmembrane GTPase MFN1
Accession No.	Swiss-Prot#: Q8IWA4NCBI Gene ID:55669
Uniprot	Q8IWA4
GeneID	55669;
Calculated MW	84 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



Western blot analysis of lysates from HCT116 cells, primary antibody was diluted at 1:1000, 4° over night

---

Note: This product is for in vitro research use only