

IBA1 Antibody

Catalog No: #45064



Package Size: #45064-1 50ul #45064-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IBA1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Rt
Specificity	IBA1 Antibody detects endogenous levels of total IBA1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human IBA1
Target Name	IBA1
Other Names	AIF1, AIF-1, Em:AF129756.17, G1, IRT-1, Protein G1, IBA1, IRT1
Accession No.	Swiss-Prot#: P55008NCBI Gene ID:199
Uniprot	P55008
GeneID	199;
Calculated MW	17 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only