

Elastin Antibody

Catalog No: #45078

Package Size: #45078-1 50ul #45078-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Elastin Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	Elastin Antibody detects endogenous levels of total Elastin
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human Elastin
Target Name	Elastin
Other Names	ELN, SVAS, WBS, Elastin, Tropoelastin, WS
Accession No.	Swiss-Prot#: P15502NCBI Gene ID:2006
Uniprot	P15502
GeneID	2006;
Calculated MW	68 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis Elastin using HepG2 whole cell lysates

Note: This product is for in vitro research use only