

## ADFP Antibody

Catalog No: #45081

Package Size: #45081-1 50ul #45081-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	ADFP Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	ADFP Antibody detects endogenous levels of total ADFP
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ADFP
Target Name	ADFP
Other Names	PLIN2, Adipophilin, ADFP, ADRP, Perilipin 2, Perilipin-2
Accession No.	Swiss-Prot#: Q99541NCBI Gene ID:123
Uniprot	Q99541
GeneID	123;
Calculated MW	48 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



Western Blot analysis of various cell lysis. Primary Antibody was diluted at 1:1000. Secondary antibody was diluted at 1:10000

---

Note: This product is for in vitro research use only