

caspase-4 Antibody

Catalog No: #45084

Package Size: #45084-1 50ul #45084-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	caspase-4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	caspase-14 Antibody detects endogenous levels of total caspase-11
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human caspase-11
Target Name	caspase-4
Other Names	CASP4; ICH2; Caspase-4; CASP-4; ICE(rel)-II; Protease ICH-2;
Accession No.	Swiss-Prot#: P49662NCBI Gene ID:837
Uniprot	P49662
GeneID	837;
Calculated MW	43 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C


Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of lysates from 1) HeLa ,2)THP-1, (Green) primary antibody was diluted at 1:1000, 4° over night, secondary antibody(cat:RS23920)was diluted at 1:10000, 37° 1hour.



Immunohistochemistry analysis of Caspase-4 antibody in paraffin-embedded human lung carcinoma tissue.

Note: This product is for in vitro research use only