

## RAP1GDS1 Conjugated Antibody

Catalog No: #C34960



Package Size: #C34960-AF350 100ul #C34960-AF405 100ul #C34960-AF488 100ul  
 #C34960-AF555 100ul #C34960-AF594 100ul #C34960-AF647 100ul  
 #C34960-AF680 100ul #C34960-AF750 100ul #C34960-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |  |
|-----------------------|--|
| Product Name          | RAP1GDS1 Conjugated Antibody   |
| Host Species          | Rabbit   |
| Clonality             | Polyclonal   |
| Species Reactivity    | Hu Ms  |
| Specificity           | The antibody detects endogenous levels of total RAP1GDS1protein.   |
| Immunogen Description | Synthesized peptide derived from C-terminal of human RAP1GDS1.   |
| Conjugates            | Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750   |
| Other Names           | KIAA0474;RAP1 GTPase activating protein;Rap1 GTPase-activating protein 1;RAP1;GTPase activating protein 1  |
| Accession No.         | Swiss-Prot#:P47736NCBI Gene ID:5909  |
| Uniprot               | P47736   |
| GeneID                | 5909;  |
| Excitation Emission   | AF350: 346nm/442nm<br>AF405: 401nm/421nm<br>AF488: 493nm/519nm<br>AF555: 555nm/565nm<br>AF594: 591nm/614nm<br>AF647: 651nm/667nm<br>AF680: 679nm/702nm<br>AF750: 749nm/775nm |
| Calculated MW         | 73   |
| Formulation           | 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide  |
| Storage               | Store at 4°C in dark for 6 months  |

## Application Details

## Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250  
 AF405 conjugated: most applications: 1: 50 - 1: 250  
 AF488 conjugated: most applications: 1: 50 - 1: 250  
 AF555 conjugated: most applications: 1: 50 - 1: 250  
 AF594 conjugated: most applications: 1: 50 - 1: 250  
 AF647 conjugated: most applications: 1: 50 - 1: 250  
 AF680 conjugated: most applications: 1: 50 - 1: 250  
 AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

---

## Product Description

---

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

## Background

---

GTPase activator for the nuclear Ras-related regulatory protein RAP-1A (KREV-1), converting it to the putatively inactive GDP-bound state.

---

Note: This product is for in vitro research use only