

SLC6A8 Conjugated Antibody

Catalog No: #C35050



Package Size: #C35050-AF350 100ul #C35050-AF405 100ul #C35050-AF488 100ul
 #C35050-AF555 100ul #C35050-AF594 100ul #C35050-AF647 100ul
 #C35050-AF680 100ul #C35050-AF750 100ul #C35050-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	SLC6A8 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total SLC6A8 protein.
Immunogen Description	Synthesized peptide derived from C-terminal of human SLC6A8.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Sodium- and chloride-dependent creatine transporter 1;CT1;Creatine transporter 1;Solute carrier family 6 member 8;SLC6A8
Accession No.	Swiss-Prot#:P48029NCBI Gene ID:6535
Uniprot	P48029
GeneID	6535;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	71
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250
 AF405 conjugated: most applications: 1: 50 - 1: 250
 AF488 conjugated: most applications: 1: 50 - 1: 250
 AF555 conjugated: most applications: 1: 50 - 1: 250
 AF594 conjugated: most applications: 1: 50 - 1: 250
 AF647 conjugated: most applications: 1: 50 - 1: 250
 AF680 conjugated: most applications: 1: 50 - 1: 250
 AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Product Description

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Background

Required for the uptake of creatine in muscles and brain.

Note: This product is for in vitro research use only