

GTF3A Conjugated Antibody

Catalog No: #C35104



Package Size: #C35104-AF350 100ul #C35104-AF405 100ul #C35104-AF488 100ul
 #C35104-AF555 100ul #C35104-AF594 100ul #C35104-AF647 100ul
 #C35104-AF680 100ul #C35104-AF750 100ul #C35104-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	GTF3A Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total GTF3A protein.
Immunogen Description	Synthesized peptide derived from Internal of human GTF3A.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	AP2;general transcription factor IIIA;GTF3A;TF3A;TFIIIA
Accession No.	Swiss-Prot#:Q92664NCBI Gene ID:2971
Uniprot	Q92664
GeneID	2971;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	42
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Product Description

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Background

Interacts with the internal control region (ICR) of approximately 50 bases within the 5S RNA genes, is required for correct transcription of these genes by RNA polymerase III. Also binds the transcribed 5S RNA's. May initiate transcription of the 5S ribosomal RNA gene and maintain the stability of transcription of other genes.

Note: This product is for in vitro research use only