

## IRG1 Antibody

Catalog No: #45116

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	IRG1 Antibody
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	IRG1 Antibody detects endogenous levels of total IRG1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human IRG1
Target Name	IRG1
Other Names	IRG1 Antibody
Accession No.	Swiss-Prot#: A6NK06NCBI Gene ID:730249
Uniprot	A6NK06
GeneID	730249;
Calculated MW	53 KD
Concentration	1.0mg/mL
Formulation	pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
Storage	Store at -20°C

## Application Details

WB 1:1000-1:5000

ELISA 1oO 5000-1oO 10000

## Images



Gel: 8%SDS-PAGE

Lysate: 40 ug

Lane: K562 cell lysate

Primary antibody: at dilution 1/1000

Secondary antibody: at 1/5000 dilution

Exposure time: 3 minutes

Note: This product is for in vitro research use only