

ITGB6 Conjugated Antibody

Catalog No: #C35783



Package Size: #C35783-AF350 100ul #C35783-AF405 100ul #C35783-AF488 100ul
 #C35783-AF555 100ul #C35783-AF594 100ul #C35783-AF647 100ul
 #C35783-AF680 100ul #C35783-AF750 100ul #C35783-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	ITGB6 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ITGB6 protein.
Immunogen Description	Fusion protein corresponding to a region derived from internal residues of human integrin, beta 6
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ITGB6, ITB6, beta 6, integrin
Accession No.	Swiss-Prot#:P18564NCBI Gene ID:3694NCBI Protein#:BC121178
Uniprot	P18564
GeneID	3694;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Background

Integrin beta-6 is a protein that in humans is encoded by the ITGB6 gene. Integrin, beta 6 has been shown to interact with FHL2. Integrin alpha-V/beta-6 is a receptor for fibronectin and cytotactin. It recognizes the sequence R-G-D in its ligands. Internalisation of integrin alpha-V/beta-6 via clathrin-mediated endocytosis promotes carcinoma cell invasion.

Note: This product is for in vitro research use only