

## PLD1 Antibody

Catalog No: #45148

Package Size: #45148-1 50ul #45148-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |   |
|-----------------------|---|
| Product Name          | PLD1 Antibody   |
| Host Species          | Rabbit  |
| Clonality             | Polyclonal  |
| Purification          | Immunogen affinity purified   |
| Applications          | WB  |
| Species Reactivity    | Hu Ms Rt  |
| Specificity           | PLD1 Antibody detects endogenous levels of total PLD1   |
| Immunogen Type        | Peptide   |
| Immunogen Description | A synthesized peptide derived from human PLD1   |
| Target Name           | PLD1  |
| Other Names           | PLD1, HPLD1, Phospholipase D1, PLD 1, Choline phosphatase 1   |
| Accession No.         | Swiss-Prot#: Q13393NCBI Gene ID:5337  |
| Uniprot               | Q13393  |
| GeneID                | 5337;   |
| Calculated MW         | 120 kDa   |
| Concentration         | 1.0mg/mL  |
| Formulation           | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide. |
| Storage               | Store at -20°C  |

## Application Details

WB (1/500 - 1/2000)

## Images



Western blot analysis of Phospholipase D1 expression in HT29 (A) whole cell lysates.

---

Note: This product is for in vitro research use only