

## GSTK1 Antibody

Catalog No: #45161

Package Size: #45161-1 50ul #45161-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	GSTK1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	GSTK1 Antibody detects endogenous levels of total GSTK1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human GSTK1
Target Name	GSTK1
Other Names	GSTK1, Glutathione S-transferase k1, HGSTK1, GST13-13, GST, GST class-kappa, GST13, GST 13-13, GSTK1-1
Accession No.	Swiss-Prot#: Q9Y2Q3NCBI Gene ID:373156
Uniprot	Q9Y2Q3
GeneID	373156;
Calculated MW	25 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



Western blot analysis of lysates from 293T cells, primary antibody was diluted at 1:1000, 4° over night

---

Note: This product is for in vitro research use only