

SCARF1 Antibody

Catalog No: #45180

Package Size: #45180-1 50ul #45180-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SCARF1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	SCARF1 Antibody detects endogenous levels of total SCARF1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human SCARF1
Target Name	SCARF1
Other Names	SCARF1, Acetyl LDL receptor, KIAA0149, SREC-I, SREC1, SREC
Accession No.	Swiss-Prot#: Q14162NCBI Gene ID:8578
Uniprot	Q14162
GeneID	8578;
Calculated MW	87 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of HeLa lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only