

## PCYT1A Antibody

Catalog No: #45181

Package Size: #45181-1 50ul #45181-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	PCYT1A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB IHC
Species Reactivity	Hu
Specificity	PCYT1A Antibody detects endogenous levels of total PCYT1A
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human PCYT1A
Target Name	PCYT1A
Other Names	PCYT1A, CCT Alpha, CCTA, CT, CTA, CCT-alpha, CT A, PCYT1, CCT A, CTPCT
Accession No.	Swiss-Prot#: P49585NCBI Gene ID:5130
Uniprot	P49585
GeneID	5130;
Calculated MW	42 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

Immunohistochemistry: 1:50~1:100

## Images



Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4° over night

---

Note: This product is for in vitro research use only