

HLA-B Antibody

Catalog No: #45189



Package Size: #45189-1 50ul #45189-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HLA-B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	HLA-B Antibody detects endogenous levels of total HLA-B
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human HLA-B
Target Name	HLA-B
Other Names	HLA-B, HLAB, MHC Class I HLA heavy chain, Leukocyte antigen class I-B
Accession No.	Swiss-Prot#: P01889NCBI Gene ID:3106
Uniprot	P01889
GeneID	3106;
Calculated MW	41 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only