

MXD3 Antibody

Catalog No: #45194



Package Size: #45194-1 50ul #45194-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MXD3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	MXD3 Antibody detects endogenous levels of total MXD3
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MXD3
Target Name	MXD3
Other Names	MXD3, BHLHC13, MAD3, MAX dimerization protein 3, MYX, Max dimerizer 3, Max-associated protein 3
Accession No.	Swiss-Prot#: Q9BW11NCBI Gene ID:83463
Uniprot	Q9BW11
GeneID	83463;
Calculated MW	23 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only