

## APOC3 Antibody

Catalog No: #45203

Package Size: #45203-1 50ul #45203-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	APOC3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	APOC3 Antibody detects endogenous levels of total APOC3
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human APOC3
Target Name	APOC3
Other Names	APOC3, Apo-CIII, Apolipoprotein C3, APOCIII, HALP2, ApoC-III, Apolipoprotein C-III
Accession No.	Swiss-Prot#: P02656NCBI Gene ID:345
Uniprot	P02656
GeneID	345;
Calculated MW	11 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1:500-2000

## Images



Western Blot analysis of 293T-UV cells using APOC3 Antibody. Primary antibody was diluted at 1:500. Secondary antibody was diluted at 1:20000.

---

Note: This product is for in vitro research use only