

RSF1 Antibody

Catalog No: #45219



Package Size: #45219-1 50ul #45219-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RSF1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	RSF1 Antibody detects endogenous levels of total RSF1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human RSF1
Target Name	RSF1
Other Names	RSF1, HBV pX associated protein-8, HBXAP, RSF-1, XAP8, HBV pX-associated protein 8, p325
Accession No.	Swiss-Prot#: Q96T23NCBI Gene ID:51773
Uniprot	Q96T23
GeneID	51773;
Calculated MW	164 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only