

OSR1 Antibody

Catalog No: #45236



Package Size: #45236-1 50ul #45236-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	OSR1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	OSR1 Antibody detects endogenous levels of total OSR1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human OSR1
Target Name	OSR1
Other Names	OSR1, ODD, Protein odd-skipped-related 1, Odd-skipped homolog
Accession No.	Swiss-Prot#: Q8TAX0NCBI Gene ID:130497
Uniprot	Q8TAX0
GeneID	130497;
Calculated MW	30 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only