

MYO18B Conjugated Antibody

Catalog No: #C37748



Package Size: #C37748-AF350 100ul #C37748-AF405 100ul #C37748-AF488 100ul
 #C37748-AF555 100ul #C37748-AF594 100ul #C37748-AF647 100ul
 #C37748-AF680 100ul #C37748-AF750 100ul #C37748-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	MYO18B Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total MYO18B protein.
Immunogen Description	Synthetic peptide of human MYO18B
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	DKFZp434C2422; DKFZp779C1668; FLJ41942; MY18B; MYO18B; myosin 18B; myosin XVIIIb; Myosin-XVIIIb
Accession No.	Swiss-Prot#:P02741NCBI Gene ID:84700NCBI Protein#:NP_000251
Uniprot	P02741
GeneID	1401;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Background

The protein encoded by this gene may regulate muscle-specific genes when in the nucleus and may influence intracellular trafficking when in the cytoplasm. The encoded protein functions as a homodimer and may interact with F actin. Mutations in this gene are associated with lung cancer.

Note: This product is for in vitro research use only