

PLA2G1B Antibody

Catalog No: #45240

Package Size: #45240-1 50ul #45240-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PLA2G1B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	PLA2G1B Antibody detects endogenous levels of total PLA2G1B
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human PLA2G1B
Target Name	PLA2G1B
Other Names	PLA2G1B, Group IB phospholipase A2, Phospholipase A2, PLA2A, PPLA2, PLA2
Accession No.	Swiss-Prot#: P04054NCBI Gene ID:5319
Uniprot	P04054
GeneID	5319;
Calculated MW	16 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of lysates from HT29, HepG2, JK, MGC 803, Mouse thymus, Mouse colon, Mouse liver, Mouse lung, Rat lung, using PLA2G1B Antibody.

Note: This product is for in vitro research use only