

## PAM Antibody

Catalog No: #45245



Package Size: #45245-1 50ul #45245-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PAM Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	PAM Antibody detects endogenous levels of total PAM
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human PAM
Target Name	PAM
Other Names	PAM, PAL, Peptidylglycine 2-hydroxylase, Peptidylamidoglycolate lyase, PHM
Accession No.	Swiss-Prot#: P19021NCBI Gene ID:5066
Uniprot	P19021
GeneID	5066;
Calculated MW	108 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only