

## SERPINA11 Conjugated Antibody

Catalog No: #C37915



Package Size: #C37915-AF350 100ul #C37915-AF405 100ul #C37915-AF488 100ul  
 #C37915-AF555 100ul #C37915-AF594 100ul #C37915-AF647 100ul  
 #C37915-AF680 100ul #C37915-AF750 100ul #C37915-Biotin 100ul

Orders: order@signalwayantibody.com  
 Support: tech@signalwayantibody.com

## Description

|                       |  |
|-----------------------|--|
| Product Name          | SERPINA11 Conjugated Antibody  |
| Host Species          | Rabbit   |
| Clonality             | Polyclonal   |
| Species Reactivity    | Hu   |
| Specificity           | The antibody detects endogenous levels of total SERPINA11 protein.   |
| Immunogen Description | Synthetic peptide corresponding to a region derived from internal residues of human serpin peptidase inhibitor, clade A (alpha-1 antiproteinase, antitrypsin), member 11     |
| Conjugates            | Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750   |
| Other Names           | ntiproteinase-like 2; Q86U17; SERPINA11  |
| Accession No.         | Swiss-Prot#:Q86U17NCBI Gene ID:256394NCBI Protein#:NP_001244960/A6NH21   |
| Uniprot               | Q86U17   |
| GeneID                | 256394;  |
| Excitation Emission   | AF350: 346nm/442nm<br>AF405: 401nm/421nm<br>AF488: 493nm/519nm<br>AF555: 555nm/565nm<br>AF594: 591nm/614nm<br>AF647: 651nm/667nm<br>AF680: 679nm/702nm<br>AF750: 749nm/775nm |
| Formulation           | 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide  |
| Storage               | Store at 4°C in dark for 6 months  |

## Application Details

## Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

## Background

---

The protein belongs to the serpin family. SERPINA11 has serine-type endopeptidase inhibitor activity, negative regulation of endopeptidase activity, regulation of proteolysis, serine-type endopeptidase inhibitor activity.

---

Note: This product is for in vitro research use only