

NEK6 Antibody

Catalog No: #45255

Package Size: #45255-1 50ul #45255-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NEK6 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	NEK6 Antibody detects endogenous levels of total NEK6
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human NEK6
Target Name	NEK6
Other Names	NEK6, NIMA-related kinase 6, Protein kinase SID6-1512, NimA-related protein kinase 6, SID6-1512
Accession No.	Swiss-Prot#: Q9HC98NCBI Gene ID:10783
Uniprot	Q9HC98
GeneID	10783;
Calculated MW	36 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only