

GAPDH Rabbit Polyclonal (Zebrafish Specific) Conjugated Antibody

Catalog No: #C38057

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #C38057-AF350 100ul #C38057-AF405 100ul #C38057-AF488 100ul

#C38057-AF555 100ul #C38057-AF594 100ul #C38057-AF647 100ul

#C38057-AF680 100ul #C38057-AF750 100ul #C38057-Biotin 100ul

Description

Product Name	GAPDH Rabbit Polyclonal (Zebrafish Specific) Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Zebrafish
Specificity	The GAPDH antibody detects Zebrafish endogenous GAPDH protein.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	G3PD; GAPD; MGC88685
Accession No.	Swiss-Prot#:P04406
Uniprot	P04406
GenID	2597;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	37
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Background

Glyceraldehyde 3 phosphate dehydrogenase (GAPDH) is well known as one of the key enzymes involved in glycolysis. GAPDH is constitutively expressed in almost all tissues at high levels, therefore antibodies against GAPDH are useful as loading controls for Western Blotting. Some physiological factors, such as hypoxia and diabetes, increase GAPDH expression in certain cell types.

Note: This product is for in vitro research use only