

## SLC9A3R1 Antibody

Catalog No: #45257

Package Size: #45257-1 50ul #45257-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	SLC9A3R1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB <sup>o</sup> ELISA
Species Reactivity	Hu Ms Rt
Specificity	SLC9A3R1 Antibody detects endogenous levels of total SLC9A3R1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human SLC9A3R1
Target Name	SLC9A3R1
Other Names	SLC9A3R1, EBP50, NPHLOP2, NHERF, NHERF-1, NHERF1
Accession No.	Swiss-Prot#: O14745NCBI Gene ID:9368
Uniprot	O14745
GeneID	9368;
Calculated MW	39 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

WB 1:500-2000 ELISA 1:5000-20000

## Images



---

Note: This product is for in vitro research use only