

## MYL4 Antibody

Catalog No: #45268

Package Size: #45268-1 50ul #45268-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	MYL4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	MYL4 Antibody detects endogenous levels of total MYL4
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MYL4
Target Name	MYL4
Other Names	MYL4, AMLC, Atrial, embryonic, ALC1, Myosin light chain 4, Myosin, light chain 4, alkali, PRO1957, GT1
Accession No.	Swiss-Prot#: P12829NCBI Gene ID:4635
Uniprot	P12829
GeneID	4635;
Calculated MW	22 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



---

Note: This product is for in vitro research use only