

ZFYVE9 Antibody

Catalog No: #45273

Package Size: #45273-1 50ul #45273-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZFYVE9 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Human
Specificity	ZFYVE9 Antibody detects endogenous levels of total ZFYVE9
Immunogen Type	Peptide
Immunogen Description	Recombinant protein of human ZFYVE9.
Target Name	ZFYVE9
Other Names	ZFYVE9, HSARA, Madh-interacting protein, Receptor activation anchor, SMADIP, MADHIP, Novel serine protease, NSP, SARA
Accession No.	Swiss-Prot#: O95405NCBI Gene ID:9372
Uniprot	O95405
GeneID	9372;
Calculated MW	156 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/2000)

Images



Western blot analysis of SARA expression in HeLa (A) whole cell lysates.

Note: This product is for in vitro research use only