

DUX4 Antibody

Catalog No: #45281

Package Size: #45281-1 50ul #45281-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DUX4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	DUX4 Antibody detects endogenous levels of total DUX4
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human DUX4
Target Name	DUX4
Other Names	DUX4, Double homeobox protein 10, Double homeobox protein DUX10, Double homeobox 4, Double homeobox protein 4/10, Double homeobox protein 4, DUX10
Accession No.	Swiss-Prot#: Q9UBX2NCBI Gene ID:100288687
Uniprot	Q9UBX2
GeneID	100288687;
Calculated MW	45 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only