

MYH2 Antibody

Catalog No: #45297



Package Size: #45297-1 50ul #45297-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MYH2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	MYH2 Antibody detects endogenous levels of total MYH2
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MYH2
Target Name	MYH2
Other Names	MYH2, Fast 2a myosin heavy chain, MyHC-IIa, Myosin heavy chain 2a, MyHC-2A, Myosin-2, MYHSA2, Myosin heavy chain 2, MYH2A, MYHas8, Myosin heavy chain IIa, Type IIA myosin heavy chain
Accession No.	Swiss-Prot#: Q9UKX2NCBI Gene ID:4620
Uniprot	Q9UKX2
GeneID	4620;
Calculated MW	223 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only