

NPC1L1 Antibody

Catalog No: #45306

Package Size: #45306-1 50ul #45306-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NPC1L1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	NPC1L1 Antibody detects endogenous levels of total NPC1L1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human NPC1L1
Target Name	NPC1L1
Other Names	NPC1L1, Niemann-Pick C1-like protein 1, NPC11L1, NPC1-like 1
Accession No.	Swiss-Prot#: Q9UHC9NCBI Gene ID:29881
Uniprot	Q9UHC9
GeneID	29881;
Calculated MW	149 kDa
Concentration	1.0mg/mL
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/2000)

Images



Western blot analysis of NPC1L1 expression in HEK293T (A), K562 (B) whole cell lysates.

Note: This product is for in vitro research use only