

## Ameloblastin Antibody

Catalog No: #45315

Package Size: #45315-1 50ul #45315-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Ameloblastin Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	Ameloblastin Antibody detects endogenous levels of total Ameloblastin
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human Ameloblastin
Target Name	Ameloblastin
Other Names	AMBN, Ameloblastin
Accession No.	Swiss-Prot#: Q9NP70NCBI Gene ID:258
Uniprot	Q9NP70
GeneID	258;
Calculated MW	48 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting:

## Images



Western blot analysis of lysates from HeLa cells, using Ameloblastin antibody.

---

Note: This product is for in vitro research use only