

RFFL Conjugated Antibody

Catalog No: #C38961



Package Size: #C38961-AF350 100ul #C38961-AF405 100ul #C38961-AF488 100ul
 #C38961-AF555 100ul #C38961-AF594 100ul #C38961-AF647 100ul
 #C38961-AF680 100ul #C38961-AF750 100ul #C38961-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | RFFL Conjugated Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Species Reactivity | Hu |
| Specificity | The antibody detects endogenous level of total RFFL antibody. |
| Immunogen Description | Recombinant protein of human RFFL. |
| Conjugates | Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750 |
| Other Names | CARP2; FRING; CARP-2; RNF189; RNF34L; RIFIFYLIN; |
| Accession No. | Swiss-Prot#:Q8WZ73NCBI Gene ID:117584 |
| Uniprot | Q8WZ73 |
| GeneID | 117584; |
| Excitation Emission | AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm |
| Calculated MW | 40 |
| Formulation | 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide |
| Storage | Store at 4°C in dark for 6 months |

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Background

#REF!

Note: This product is for in vitro research use only